

Defibrillator checklist

In the UK, the NHS ambulance services attempt resuscitation in approximately 30,000 people each year. Bystander CPR and use of a defibrillator can increase the chances of survival by two to fourfold*.

We've put together a checklist to help you ensure your defibrillator is always ready for use.

Defibrillators self-test on a daily, weekly or monthly basis - so your device will signal, typically with a flashing light or audible alert, if there is a problem, for example, pads not correctly connected, or low battery.

1	

Set up a regular inspection procedure for your defibrillator



Record your inspections below




Record the expiration date of pads and **batteries**



Reorder batteries and pads before they expire



Follow the guidance in your product manual when cleaning your defibrillator



Defibrillator ready to use?	Date	Initial



Need new pads or batteries? Order now at sja.org.uk/shop or call 0344 770 4808†

* 2021 Resuscitation Council Guidelines, Resuscitation Council UK